

Mini-spectrometer TM series

C10082MD, C10083MD

Integrating optical system, image sensor and circuit



TM series mini-spectrometers are polychromators integrated with optical elements, an image sensor and a driver circuit. Two models are available: C10082MD (TM-UV/VIS-MOS) and C10083MD (TM-VIS/NIR-MOS). Light to be measured is guided into the entrance port of TM series through an optical fiber and the spectrum measured with the built-in image sensor is output from the USB port to a PC for data acquisition. No external power supply is required since USB bus power is used for circuit operation. The TM series comes with sample software, device driver and DLL. The sample software lets you easily set measurement conditions, acquire and save data, and display graphs.

Features

- High throughput due to transmission grating made of quartz
- Highly accurate optical characteristics
- No external power supply required: Uses USB bus power
- Wide spectral response range
- Easy to install into equipment
- Wavelength conversion factor *1 is recorded in internal memory

Applications

- Light source spectrum measurement
- Sunlight or illumination analysis
- Absorbance measurement

■ Optical characteristics

| Parameter | TM-UV/VIS-MOS | TM-VIS/NIR-MOS | Unit |
|---|---------------|----------------|-------|
| | C10082MD | C10083MD | |
| Spectral response range | 200 to 800 | 320 to 1000 | nm |
| Spectral resolution Max. (Spectral response half width) *2 | 6 | 8 | nm |
| Wavelength reproducibility *3 | | ±0.2 | nm |
| Wavelength temperature dependence | | 0.02 | nm/°C |
| Spectral stray light *2, *4 | | -35 | dB |

■ Electrical characteristics

| Parameter | Value | Unit |
|-----------------------------------|------------|------|
| A/D conversion | 16 | bit |
| Integration time | 5 to 10000 | ms |
| Interface | USB1.1 | - |
| USB bus power current consumption | 100 | mA |

■ General ratings / Absolute maximum ratings

| Parameter | Value | Unit |
|-----------------------------|--|--------|
| Dimensions | 94 (W) × 90 (D) × 55 (H) | mm |
| Image sensor | CMOS linear image sensor (S8378-1024Q) | - |
| Number of pixels | 1024 | pixels |
| Slit *5 | 70 (H) × 800 (V) | μm |
| Optical NA | 0.22 | - |
| Connector for optical fiber | SMA905D | - |
| Operating temperature *6 | +5 to +40 | °C |
| Storage temperature | -20 to +70 | °C |

*1: A conversion factor for converting the image sensor pixel number into a wavelength is recorded in the module. A calculation factor for converting the A/D converted count into the input light intensity is not provided.

*2: Depends on the slit opening. Values were measured with the slit opening listed in the table.

*3: Measured under constant light input conditions

*4: When monochromatic light of the following wavelengths is input, spectral stray light is defined as the ratio of the count measured at the input wavelength, to the count measured in a region of the input wavelength ±40 nm.
C10082MD: 500 nm, C10083MD: 650 nm

*5: Entrance slit aperture size

*6: No condensation


Light spectrally separated by a grating is focused according to light wavelength on predetermined image sensor positions, and high-order light is also input onto positions at integer-multiples of wavelengths. In these mini-spectrometers an optical filter is attached to the image sensor to cut off high-order light, but this also causes a drop in the image sensor output at the following wavelengths.

C10082MD: Near 340 nm and 500 nm, C10083MD: Near 500 nm and 700 nm

Types not using a high-order light cut off filter are also available. Please specify by adding "-01" to the type number when ordering. (Example: C10082MD_01)

Optical component layout


TM series mini-spectrometers use a transmission holographic grating made of quartz and precision optical components arranged on a rugged optical base, making it possible to deliver high throughput and highly accurate optical characteristics.



KACCC0287EA

Connection example (transmission light measurement)

Light to be measured is guided into the entrance port of TG series through an optical fiber and the spectrum measured with the built-in image sensor is output through the USB port to a PC for data acquisition. There are no moving parts inside the unit so stable measurements are obtained at all times. An optical fiber that guides light input from external sources allows a flexible measurement setup.



KACCC0288EF


Dedicated software package (supplied with unit)

Installing the dedicated software package (containing sample software, device driver, DLL)*7 into your PC allows running the following basic tasks:

- Measurement data acquisition and save
- Measurement condition setup
- Module information acquisition (wavelength conversion factor, polychromator type, etc.)
- Graphic display
- Arithmetic operation
 - Pixel number to wavelength conversion
 - Dark subtraction
 - Comparison calculation with reference data (transmittance, reflectance)
 - Gaussian approximation (peak position and count, FWHM)

Note: Two or more mini-spectrometers can be connected and used with one PC simultaneously.

*7: Compatible OS: Microsoft Windows Professional Edition 2000 (SP3 or later) and XP (SP1a or later)



Device driver and DLL for controlling hardware are also provided.

You can develop your own measurement programs by using a software development environment that includes Microsoft Visual C++ and Visual Basic.*8 The DLL provides functions such as USB port open/close, measurement condition setup, measurement data and module information acquisition.


*8: Operation of the device driver and DLL has been verified only with Microsoft Visual C++® and Visual Basic®.

Microsoft Visual C++ and Microsoft Visual Basic are either registered trademarks or trademarks of Microsoft Corporation in the United States.


Mini-spectrometer TM series C10082MD, C10083MD

■ Measurement example

Line spectra from LED were measured with C10082MD (TM-UV/VIS-MOS)



■ Dimensional outline (unit: mm)



KACCB0126EA

KACCA0171EC

■ Accessories

- USB cable
- Dedicated software (sample software, device driver, DLL)

■ Options (sold separately)

Optical fibers for light input

| Type. No | Product name | Core diameter (μm) | Specification |
|----------|---|--------------------|---|
| A9762-01 | Fiber for UV/visible range (resistance to UV) | 600 | N.A.=0.22, length 1.5 m, connectorized SMA905D at both ends |

Mini-spectrometer TM series C10082MD, C10083MD

■ Mini-spectrometer line-up

| Type No. | Type | Spectral response range (nm) | | | | | | | | | | | Spectral resolution Max. (nm) | Image sensor | | | | |
|-----------|---|------------------------------|-------------|-----|-------------|------|------|------|------|------|------|------|-------------------------------|--------------|--|-------------------------|---------------------------------------|-----------------------------|
| | | 200 | 400 | 600 | 800 | 1000 | 1200 | 1400 | 1600 | 1800 | 2000 | 2200 | | | | | | |
| C10082CA | TM-UV/VIS-CCD High sensitivity | | | | | | | | | | | | | | | 6 | Back-thinned type CCD image sensor | |
| C10082CAH | TM-UV/VIS-CCD High resolution | | 200 to 800 | | | | | | | | | | | | | 1* | | |
| C10082MD | TM-UV/VIS-MOS Wide dynamic range | | | | | | | | | | | | | | | 6 | CMOS linear image sensor | |
| C10083CA | TM-VIS/NIR-CCD High sensitivity | | | | | | | | | | | | | | | 8 (λ=320 to 900 nm) | Back-thinned type CCD image sensor | |
| C10083CAH | TM-VIS/NIR-CCD High resolution | | 320 to 1000 | | | | | | | | | | | | | 1* (λ=320 to 900 nm) | | |
| C10083MD | TM-VIS/NIR-MOS Wide dynamic range | | | | | | | | | | | | | | | 8 | CMOS linear image sensor | |
| C9404CA | TG-UV-CCD High sensitivity | | | | | | | | | | | | | | | 3 | Back-thinned type CCD image sensor | |
| C9404CAH | TG-UV-CCD High resolution | | 200 to 400 | | | | | | | | | | | | | 1* | Back-thinned type CCD image sensor | |
| C9404MC | TG-UV-MOS Wide dynamic range | | | | | | | | | | | | | | | 3 | CMOS linear image sensor | |
| C9405CA | TG-SWNIR-CCD High sensitivity | | | | | | | | | | | | | | | 5 (λ=550 to 900 nm) | Back-thinned type CCD image sensor | |
| C9405MC | TG-SWNIR-MOS Wide dynamic range | | | | 500 to 1100 | | | | | | | | | | | 5 (λ=550 to 1100 nm) | | |
| C9406GC | TG-NIR Non-cooled type | | | | | | | | | | | | | | | 7 | InGaAs linear image sensor | |
| C9913GC | TG-cooled NIR-I Low noise (cooled type) | | | | | | | | | | | | | | | 7 | | |
| C9914GB | TG-cooled NIR-II Low noise (cooled type) | | | | | | | | | | | | | | | 8 | | |
| C9407MA | RC-VIS-MOS Spectrometer module | | 340 to 780 | | | | | | | | | | | | | | 9 | CMOS linear image sensor |

* Typ.

OEM model

| Type No. | Type | Spectral response range (nm) | | | | | | | | | | | Spectral resolution Max. (nm) | Image sensor | | | | |
|----------|---------------------------------|------------------------------|------------|-----|-----|------|------|------|------|------|------|------|-------------------------------|--------------|--|--|---|-----------------------------|
| | | 200 | 400 | 600 | 800 | 1000 | 1200 | 1400 | 1600 | 1800 | 2000 | 2200 | | | | | | |
| C9409MA | RC-VIS-MOS Spectrometer head | | 340 to 780 | | | | | | | | | | | | | | 9 | CMOS linear image sensor |

HAMAMATSU

Information furnished by HAMAMATSU is believed to be reliable. However, no responsibility is assumed for possible inaccuracies or omissions. Specifications are subject to change without notice. No patent rights are granted to any of the circuits described herein. ©2007 Hamamatsu Photonics K.K.

HAMAMATSU PHOTONICS K.K., Solid State Division

1126-1 Ichino-cho, Higashi-ku, Hamamatsu City, 435-8558 Japan, Telephone: (81) 53-434-3311, Fax: (81) 53-434-5184, www.hamamatsu.com

U.S.A.: Hamamatsu Corporation: 360 Foothill Road, P.O.Box 6910, Bridgewater, N.J. 08807-0910, U.S.A., Telephone: (1) 908-231-0960, Fax: (1) 908-231-1218

Germany: Hamamatsu Photonics Deutschland GmbH: Arzbergerstr. 10, D-82211 Herrsching am Ammersee, Germany, Telephone: (49) 08152-3750, Fax: (49) 08152-2658

France: Hamamatsu Photonics France S.A.R.L.: 19, Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: 33-(1) 69 53 71 00, Fax: 33-(1) 69 53 71 10

United Kingdom: Hamamatsu Photonics UK Limited: 2 Howard Court, 10 Tewin Road, Welwyn Garden City, Hertfordshire AL7 1BW, United Kingdom, Telephone: (44) 1707-294888, Fax: (44) 1707-325777

North Europe: Hamamatsu Photonics Norden AB: Smidesvägen 12, SE-171 41 Solna, Sweden, Telephone: (46) 8-509-031-00, Fax: (46) 8-509-031-01

Italy: Hamamatsu Photonics Italia S.R.L.: Strada della Moia, 1/E, 20020 Arese, (Milano), Italy, Telephone: (39) 02-935-81-733, Fax: (39) 02-935-81-741

Cat. No. KACC1119E05
Jan. 2007 DN